

Finnigan LTQ FT

Ultra High Performance Mass Spectrometer

Finnigan LTQ FT – Ultra High Performance

- Ion Trap based Fourier Transform ICR Mass Spectrometer
- FTICR (FTMS) at LC timescale
- Accurate Mass
- High Sensitivity
- Ultra-high Resolution
- Above All - **Simple to Operate**



Finnigan LTQ FT – Ultra High Performance

- Surveyor HPLC
- API, MALDI Sources
- 2-D Linear Ion Trap
- ICR Cell
- Actively Shielded Superconducting Magnet
- **For details – have a look at the dummy**




Finnigan LTQ FT – Ultra High Performance

- First Routine FTMS Mass Spectrometer
 - *Fast*
 - *Automated*
 - *Data dependent scanning*
 - *Robust technology*
- Most Advanced Mass Spectrometer Available
 - *Proteomics*
 - *Metabolomics*
 - *Pharmaceutical and Drug Discovery*
 - *Small Molecule Analysis, Structural Elucidation*

Finnigan LTQ FT – Combination of Performance

- All Useful Ion Trap Features
 - API with APCI, ESI and MALDI Ion Sources
 - MS, MS/MS and MSⁿ Analysis
 - AGC Control
 - Data Dependent Scanning
 - Secondary Electron Multiplier *Detector*



Finnigan LTQ FT – Experiment Description

- **Generate Ions**
 - *API: ESI, APCI, APPI, MALDI Source*
- **Store and Prepare Ions in 2-D Linear Trap**
 - *MS, MSⁿ*
 - *AGC to regulate number of ions*
- **Transfer Ions to ICR Cell**
 - *RF ion guides*
- **Excite Ions in ICR Cell**
 - *RF frequency sweep or SWIFT excitation*
- **Detect Ion Transient Signal**
- **Fourier Transform Data Processing**
 - *Data dependent information passed to 2-D Linear Trap*

Finnigan LTQ FT – Performance Specifications

- **Resolution**
 - *100 000 resolution at m/z 400 at 1 Hz repetition rate*
 - *>500 000 resolution broadband mode*
- **Mass Range**
 - *m/z 50-4000 (standard range)*
 - *1-order-magnitude in single scan (e.g. m/z 200-2000)*
- **Mass Accuracy**
 - *2 ppm RMS, external mass calibration*
 - *<1 ppm RMS, internal mass calibration*
- **Dynamic Range**
 - *5 000 within mass spectrum*
 - *>500 000 between mass spectra*

Finnigan LTQ FT – Accurate Mass

- The magnetic field is very stable, therefore mass calibration is valid for several days
- Accurate Mass with External Calibration
 - *Must control the number of ions in ICR cell - Ion cyclotron frequency varies with ion density*
 - *Number of ions controlled with AGC on the 2-D Linear Trap*
 - *Mass accuracy of 2 ppm RMS for MS and MS/MS mode*
- Accurate Mass with Internal Calibration
 - *A reference substance of known mass is analyzed with the unknown substance and used as a mass calibrant*
 - *Mass accuracy of better than 1 ppm RMS*

Finnigan LTQ FT – High Resolution

- Resolution Varies with Mass
 - Resolution is inversely proportional to mass
- Resolution Depends on Detection Time
 - Double the detection time to double the resolution
 - Detect time of 745 ms for 100,000 resolution at m/z 400
- Resolution Depends on Magnet Field Strength
 - Double the magnet field to double the resolution
 - Magnet price increases immensely with field strength !

Finnigan LTQ FT – 100 000 Resolution at 1 Hz

Calibration Mixture
m/z 195 to 1922


1024 k Word transient
Acquisition rate 1 Hz

MRFA – ^{13}C and ^{34}S Isotopes resolved

Res. Power at m/z 526:
90 000 (FWHM)


This is equivalent to

Res. Power at at m/z 400:
118,000 (FWHM)





Finnigan LTQ FT – Ultra High Resolution

NRVYVHPFHL (3+) 427.89988 (Angiotensine from Goosefish)



Finnigan LTQ FT – MS/MS Analysis


MS/MS Angiotensin I DRVYIHPFHL m/z 432.9 (3+), single Scan



Finnigan LTQ FT – Productive Software

- **Known and Accepted Software Platform**
 - *LCQ TUNE PLUS enhanced for FTICR analysis*
 - *Data dependent scanning*
 - *Xcalibur enhanced for accurate mass*
 - *Bioworks*
 - *DeNovoX*
 - *Metabolite ID*
 - *Mass Frontier*

Xcalibur™ Software




TunePlus

- Simple to use
- Intuitive GUI

Xcalibur 1.4


- Extended functions for high resolution and accurate mass



Finnigan LTQ FT – Productive Software

Known Platform

- LCQ **TUNE PLUS**, enhanced for FTICR analysis
- Method Editor, with data dependent scanning
- **Xcalibur**
- Bioworks, DeNovoX, Metabolite ID, Mass Frontier
- Accurate Mass Data Evaluation Software



If you can run a Finnigan LCQ you can run an LTQ FT

LTQ FT – Tune Plus Window

The screenshot shows the 'Tune Plus' software interface for a Finnigan LTQ FT instrument. The main window is titled 'Define Scan' and contains several sections:

- Scan Description:** Analyzer (FTMS), Mass Range (High), Resolution (12500), Scan Type (12500), Scan Time (100000), Microscans (400000), Mag. Inject Time (ms) (50.000).
- Source Fragmentation:** Qn, Energy (V): 20.0.
- MSn Settings Table:**

| n | Parent Mass (m/z) | Isolation Width (m/z) | Normalized Collision Energy | Activation Q | Activation Time (ms) |
|----|-------------------|-----------------------|-----------------------------|--------------|----------------------|
| 2 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 3 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 4 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 5 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 6 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 7 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 8 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 9 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 10 | | 1.0 | 20.0 | 0.250 | 30.000 |

- Scan Ranges Table:**

| # | First Mass (m/z) | Last Mass (m/z) |
|----|------------------|-----------------|
| 1 | 150.00 | 2000.00 |
| 2 | | |
| 3 | | |
| 4 | | |
| 5 | | |
| 6 | | |
| 7 | | |
| 8 | | |
| 9 | | |
| 10 | | |

Buttons at the bottom of the 'Define Scan' window include 'Apply', 'OK', 'Cancel', 'Help', 'Injection RE...', and 'Activation...'. A red arrow points to the 'Define Scan Button' label above the 'Define Scan' window. Another red arrow points to the 'Cryo Button' label above the 'FT Cryo Monitor' dialog box. A third red arrow points to the 'ICR Cell Button' label above the 'Injection Control' dialog box.

Define Scan Button

FT Cryo Monitor

- Liquid Level
- Helium (%):
- Nitrogen (%):

Cryo Button

Close Help

Injection Control

Ion Trap: FT

AGC Settings:

- Full MS Target: 200000.0
- SIM Target: 100000.0
- MSn Target: 100000.0
- Dyn Target: 100000.0

Enable Injection Waveforms

OK Cancel Apply Help

Method Setup for Data Dependant Scanning

Method Editor

ft_double_1st.meth - Instrument Setup

File LTQ Help

Endeavor
LTQ MS

MS Detector Setup | Syringe Pump | Divert Valve | Contact Closure | Summary

Run settings
Acquire time (min): 0.50 Segments: 1 Start delay (min): 0.00

To display a chromatogram here, use Jupiter/Open raw file...

Segment 1
Retention time (min)

Segment 1 settings
Segment time (min): 0.50 Scan events: 2 Tune method: C:\Xcalibur\system\ltq\msx\tune_pt1_...


Scan Event 1 Scan Event 2

Scan event 2 settings
Analyzer: Ion Trap FT
Mass range: Normal High Low
Scan mode: MS MS/MS MSn Parent masses...
Scan type: Full SRM
Scan rate: Normal Turbo Zoom UltraZoom
Mass Range: Setup ranges...
Polarity: Positive Negative
Data type: Centroid Profile

Source fragmentation
 Turn on Energy (V): 20.0
 Wideband activation
 Dependent scan Settings...


New method Tune Plus Help

Ready



Additions made to standard method editor

MS/MS of MRFA – Scan Setup and Tune Plus



Additions made to standard 'Define Scan' editor

Define Scan Button

Define Scan

Scan History: [Dropdown]

Scan Description:

- Analyzer: FTMS
- Mass Range: Normal
- Resolution: 100000
- Scan Type: Full

Scan Time:

- Microscans: 1
- Max. Inject Time (ms): 500.000

Source Fragmentation:

- Qn Energy (V): 20.0
- Wideband Activation

MSn Settings:

| n | Parent Mass (m/z) | Isolation Width (m/z) | Normalized Collision Energy | Activation Q | Activation Time (ms) |
|----|-------------------|-----------------------|-----------------------------|--------------|----------------------|
| 2 | 524.00 | 5.0 | 30.0 | 0.250 | 30.000 |
| 3 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 4 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 5 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 6 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 7 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 8 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 9 | | 1.0 | 20.0 | 0.250 | 30.000 |
| 10 | | 1.0 | 20.0 | 0.250 | 30.000 |

Scan Ranges:

| # | First Mass (m/z) | Last Mass (m/z) |
|----|------------------|-----------------|
| 1 | 150.00 | 2000.00 |
| 2 | | |
| 3 | | |
| 4 | | |
| 5 | | |
| 6 | | |
| 7 | | |
| 8 | | |
| 9 | | |
| 10 | | |

Input: From/To

Buttons: Apply, OK, Cancel, Help, Injection RE..., Activation...

Finnigan LTQ FT – Comparison

- Unlike Other FTICR Instruments...
 - *no trapping gas is pulsed into the ultra-high vacuum region to trap ions in the ICR cell*
 - don't need to wait several seconds to pump away gas before detecting ICR signal
 - *no collision gas is pulsed into the ICR cell for MS/MS*
 - MS/MS performed in 2-D Linear Trap
 - again, don't need to wait to pump away gas prior to detection
 - *minimized transmission discrimination*
 - short path length between 2-D Linear Trap and ICR cell and octopole ion optics, reduces mass discrimination
 - *ultra High Speed Signal Processing*
 - Fourier Transformation performed in real time

Finnigan LTQ FT – What About ?

- High Resolution Parent Ion Isolation?
 - *possible in ICR cell with SWIFT isolation*
- Collision Activated Dissociation in the ICR Cell?
 - *need pulsed gas inlet. Not planned*
 - *MS/MS performed in 2-D Linear Trap*
- IRMPD?
 - *possible, and in progress*
- ECD?
 - *possible, and in progress*

High Resolution and Accurate Mass

Peptides: [Val⁵]-Angiotensin II
Sequence: DRVYVHPF
Formula: C₄₉H₆₉N₁₃O₁₂
Exact mass: [M+2H]²⁺ = 516.76671

Lys-des-Arg⁹-Bradykinin
KRPPGFSPF
C₅₀H₇₃N₁₃O₁₁
[M+2H]²⁺ = 516.78490

Δm (mmu):



18.2 mmu

Δm (ppm):

35 ppm

Required Resolution:

56,700 (FWHM)



Mass Error and Elemental Composition


Entering the exact mass of a peptide into a program to calculate the elemental composition will yield different numbers of proposals. These numbers only depend on the maximum allowed mass errors and limits for the elements in use. The below listed limits of elements is for an average peptide composed of common amino acids. The mass error differs with instrument type and is not directly related to the resolution if a single compound is analyzed.

Example Peptide: [Val⁵]-Angiotensin II, [M+2H]²⁺ = 516.77671


Limitation for Elements: C 35-70, H 45-100, N 8-16, O 9-16, S 0-2

| Instrument | | | QqTOF | LTQ-FT | LTQ-FT |
|----------------------------------|--|--------|-------|--------|--------|
| Mass Error | | 10 ppm | 5 ppm | 2 ppm | 1 ppm |
| | | | | | |
| # of Proposals for m/z 516.76671 | | 49 | 23 | 10 | 4 |


Finnigan LTQ FT – Angiotensin I



432.89977
(error -0.05 ppm)




Sample Concentration: 40 fmol/ μ l
Nano-spray Flow Rate: ~ 30 nl/min
Sample Consumption: ~ 20 amol/s
Single Scan
Injection Time: ~ 1 ms



648.84602
(error -0.3 ppm)


Finnigan LTQ FT – MS/MS Angiotensin I

T: + p ESI Full ms2 433.00@28.00 [195.00-2000.00]




Finnigan LTQ FT – MS/MS Angiotensin I

+ p ESI Full ms2 433.00@28.00 [195.00-2000.00]




Finnigan LTQ FT – MS³ Angiotensin I

+ p ESI Full ms3 433.00@32.00 583.00@14.00 [195.00-2000.00]



Using the two MS detectors (Linear ion trap and FTMS) for post-translational modification assignment




Using the two MS detectors to pinpoint the PTM – phosphopeptides



Using the two MS detectors to pinpoint the PTM – phosphopeptides

The neutral loss of the previous ion (1031 m/z) is the major peak (982 m/z) on MS/MS spectrum as shown. Using the accurate mass and high resolution on the FTMS in MS/MS mode (confirmed that's a doubly charge ion with the loss of H₃PO₄ – it's not the loss of Pro or Val or the loss of sialic acid in 3+ charges)




Tryptic peptide mixture



Using the two MS detectors to pinpoint the PTM – disulfide-linked peptides

Using the accurate mass and high resolution on FTMS full MS mode (0.1 ppm mass accuracy and 3+ charge) confirmed the two cysteine-containing peptides are linked

Tryptic peptide mixture




Time: 12.12, Intensity: 6.06e+003, Scan Filter: (none)

NUM

Using the two MS detectors (Linear ion trap and FTMS) for intact protein assignment

Human Growth Hormon (intact): FTMS Spectrum



Finnigan LTQ FT

Recent Results


AP-MALDI on the Finnigan LTQ FT

FPPHHLTIVLGR
1386.80047
0.4 ppm

MALDI_030524182253 # 390 RT: 8.28 AV: 1 NL: 3.62E3
T: FTMS + p ESI Full ms [200.00-2000.00]


MALDI

Full scan MS
20 Laser shots
Masstech AP-MALDI
RP= 100,000







AP-MALDI MS/MS with the Finnigan LTQ FT

MALDI_030524182253 #418-438 RT: 10.68-12.08 NL: 1.72E2
FTMS + p ESI w Full ms2 1386.80@45.00 [200.00-2000.00]



LC-MS Accurate Mass with the LTQ FT









Negative ESI



Accurate Mass for Small Molecules


On-line HPLC measurements of 4 different compounds were performed on the LTQ FT at a resolution power of 200,000. AGC was used to control the number of ions.

| Cinchonidine | | | |
|---------------------|----------|-----------|--|
| Scan 55 to 155 | | | |
| Accurate mass: | | 295.18049 | |
| Aver: | 295.1802 | -0.91 ppm | |
| Min: | 295.1801 | -1.24 ppm | |
| Max: | 295.1804 | -0.21 ppm | |
| Lidoflazine | | | |
| Scan 521 to 551 | | | |
| Accurate mass: | | 492.28210 | |
| Aver: | 492.2816 | -1.02 ppm | |
| Min: | 492.2813 | -1.74 ppm | |
| Max: | 492.2820 | -0.31 ppm | |
| Procaïnamide | | | |
| Scan 55 to 155 | | | |
| Accurate mass: | | 236.17574 | |
| Aver: | 236.1756 | -0.51 ppm | |
| Min: | 236.1755 | -1.1 ppm | |
| Max: | 236.1757 | 0.0 ppm | |
| Trimipramine | | | |
| Scan 480 to 630 | | | |
| Accurate mass: | | 295.21688 | |
| Aver: | 295.2166 | -0.97 ppm | |
| Min: | 295.2165 | -1.34 ppm | |
| Max: | 295.2168 | -0.31 ppm | |



Accurate Mass for Small Molecules

Resolution = 200,000



2 ppm window


TIC

NL: 4.93E6
m/z = 295.17990-295.18108


NL: 1.60E7
m/z = 295.21629-295.21747

NL: 2.34E6
m/z = 236.17527-236.17621


NL: 6.61E6
m/z = 492.28112-492.28308




Procainamide



Cinchonidine




Trimipramine




Lidoflazine

The Advion NanoMate 100 on the Finnigan LTQ FT




Mass Accuracy for Repetitive Measurements




Mass accuracy for repetitive measurements of the MH⁺ Peak of Bupropion on the LTQ FT with the NanoMate 100 attached . 13 injection were carried out giving the 2400 data points. Same results are obtained for 3 other small molecules, Threhydrobupropion, Hydroxybupropion (HB), D₆-BUP and D₆-HB. The red line indicates the theoretical mass. **The average error is -0.20ppm.**

Accurate Mass with the LTQ FT


The combination of the Advion Nanomate system with the LTQ FT mass spectrometer was used to make 13 repetitive measurements of Bupropion. Statistics for the molecular ion and its most abundant isotope peaks are shown in the graph. For each injection the Nanomate performed a nanospray at ca. 200 nL/min. The complete set of data was acquired in less than 1 hour.




Accurate Mass Statistics for Bupropion M+H



Accurate Mass Statistics for Bupropion M+H+1




Accurate Mass Statistics for Bupropion M+H+2




The Advion NanoMate 100 on the Finnigan LTQ FT

Variation of Sample Concentration for the Measurement of Threohydrobupropion





Threohydrobupropion
 $[M+H]^+ = 242.13062$
 $C_{13}H_{20}NOCl$



The acquired spectra were analyzed for the accurate masses of the four isotope peaks on the protonated molecule. The mass chromatogram for the MH^+ ion of Threohydrobupropion is shown as well as one exemplary spectrum. Again, the first four isotope peaks on the MH^+ ion were evaluated. As expected, the mass accuracy for all different concentrations are equally good and there is no trend visible with regard to mass accuracy. The errors are randomly distributed across the measurements, all results are better than 0.5 ppm.

Ubiquitin at Very High Resolution

Settings: 1e-6 M, infusion @ 1 ul/min, 5 transients averaged, 4.5 ms inject time / transient




Ubiquitin at Very High Resolution – A Comparison


Settings: 1e-6 M, infusion @ 1 ul/min, 5 transients averaged

| | | |
|-------------------------------------|--------------------|-------------------|
| `IonSpec Experiment`: (Web page) | Resolution | 200,000 @ m/z 779 |
| | Sample consumption | 30 fmol |
| | S/N | 115 : 1 |
| | Inject time | 360 ms |

| | | |
|------------------|--------------------|-------------------|
| Finnigan LTQ-FT: | Resolution | 200,000 @ m/z 779 |
| | Sample consumption | 408 amol |
| | S/N | 270 : 1 |
| | Inject time | 24.5 ms |






Blow up →



On-line Nano HPLC-MS at Ultra High Resolution

Sample: very complex crude cell extract (non digested) from human blood platelets
Amount: unknown, but very low conc.
Flow: 200 nl / min

Resolution:
400,000 @ m/z 400



Blow up

Workshop

Finnigan LTQ FT

Thermo Electron's Annual Customer Forum
Sunday, June 9, 2003
Montreal, Canada

Agenda

- **Welcome and Introduction**
Helmut Muenster, Thermo Electron Corp., Bremen, Germany
- **Finnigan LTQ FT – A Technical description of the new Hybrid FTMS**
Stevan Horning, Thermo Electron Corp., Bremen, Germany
- **High Resolution FTMS for Quantitative Proteomics Measurements**
Thomas P. Conrads, National Cancer Institute (NCI) Frederick, NC
- **First Results with the new LTQ FT**
Michael Linscheid, Humboldt University, Berlin, Germany



Demo Chemist – Jens Griep-Raming



Jens Griep-Raming
Bremen, Germany
+49-421-5493-219

Poster MPI 182
Fully Automated High-Throughput Accurate
Mass Determination using FT-ICR Mass
Spectrometry

Jens.Griep-Raming@Thermo.com

Demo Chemist – Wolfgang Metelmann-Strupat



Wolfgang Metelmann-Strupat
Bremen, Germany
+49-421-5493-255

Poster TPX 449

An Integrated LC-MS Platform for Rapid and Precise Identification of Proteins and their Post-Translational Modifications using a Linear Ion Trap coupled with a FT-ICR Mass Spectrometer

Wolfgang.Metelmann-Strupat@Thermo.com

Demo Chemist – Billy Wu



Shiao-Lin (Billy) Wu
San Jose, USA
+1-408-965-6308

Poster TPE 091
Fourier Transform Ion Cyclotron Resonance
(FT-ICR) Mass Spectrometry – A Rapid
Method for Metabolite Identification

Billy.Wu@Thermo.com

Additional LTQ FT Posters at ASMS

Poster MPI 189

Timo Hagemeister, Humboldt Univ. Berlin, Germany

**Fragmentation Pathways of Cisplatin Adducts to Dinucleotides
Determined by FT-ICR-MS**

Poster MPX 457

Michael Linscheid, Humboldt Univ. Berlin, Germany

**Structure Elucidation of Structural Proteins from Yersina Phages
using MALDI-ToF and ESI-FTMS Data**

Finnigan LTQ FT

Installation Requirements

Superconducting Magnet – Requirements

- **Transportation**
 - *Cold shipping possible*
- **Positioning in Laboratory**
 - *Active shielding so no influence on lab equipment*
 - *Must not be positioned near elevator, large electric motor or massive steel structural supports*
- **Maintenance**
 - *Autofill for LN2*
- **Safety**
 - *Venting*
 - *Oxygen depletion sensors*


Finnigan LTQ FT

Conclusion

Highlights

- High resolution @ 1 Hz repetition rate
- Automation
- High mass accuracy
- On-line LC-MS
- MS, MS/MS and MSⁿ
- Low maintenance
- Small footprint
- And overall – **simple to use**

Bremen Development Team



... plus R&D Team in San Jose !!!