

# The ThermoFinnigan TSQ Quantum

**High Performance Quantitative and Qualitative Applications  
on a High Resolution Triple Quadrupole Mass Spectrometer**

**Gary Paul, Witold Winnik, Maurizio Splendore**

**Thermo Finnigan**


# The TSQ Quantum

The **only** high resolution triple quadrupole mass spectrometer available capable of accurate mass measurement.

Identification and high performance quantitation of components such as metabolites in one **small** package!



# TSQ Quantum - How So Small? A Revolutionary 90 Degree Bent Collision Cell




# TSQ Quantum - Quantitation

High performance quantitation is the “bread and butter” application of the triple quadrupole mass spectrometer market

TSQ Quantum was primarily designed to achieve maximum **sensitivity**, **accuracy**, **precision** and **linear dynamic range**

# API MS/MS Selective Reaction Monitoring (SRM)




# TSQ Quantum SRM Quantitative Assays at Unit Resolution

**Sensitivity**  
**What Everyone Wants**

# TSQ Quantum, Low Femtogram +ESI/SRM at Unit Resolution (Q1, Q3 - 0.7 Da FWHM)



**Proprietary Pharmaceutical** in Mobile Phase, 3 $\mu$  Phenyl Hexyl 20 x 2mm,  
ACN/0.1%TFA/H<sub>2</sub>O/NH<sub>4</sub>CH<sub>3</sub>COO, +ESI, SRM m/z 425 - 218, Collision Energy 28V



# TSQ Quantum


## Low Femtogram +ESI/SRM at Unit Resolution

Descarboethoxyloratadine in Human Plasma, H<sub>2</sub>O/MeOH/ACN/0.1% formic,  
250  $\mu$ l/min, 100 x 2 mm C8, +ESI, SRM m/z 311 - 259, Collision Energy 28V




# TSQ Quantum, Low Femtogram +ESI/SRM at Unit Resolution

## Plasma Blank



## 10 fg on column **Formoterol**




Thanks to Mark Churchill and Mark Harrison  
for data

**Thermo**Finnigan

# TSQ Quantum – Nonapeptide (MW 1209) in Human Plasma 0.8 fM on column, +ESI/SRM in High Resolution

C18 50 mm x 4.6 mm, ACN/0.1%FA/H<sub>2</sub>O/0.1%FA, SRM m/z 605 - 249, CE 30V


SM: 7G



# TSQ Quantum SRM Quantitative Assays at Unit Resolution

**Why So Sensitive?**

# The New TSQ Quantum Analyzer



# TSQ QUANTUM - Reasons for Improved S/N

1. **New Orthogonal API Source - Less Noise**
2. **New Square Quadrupole Ion Guides – More Signal**
3. **New Hyperbolic Quadrupoles - More Signal**
4. **New 90 Degree Collision Cell –  
More Signal, Less Noise**

# TSQ Quantum

## SRM Quantitative Assays at Unit Resolution


### ALSO ESSENTIAL FOR HIGH PERFORMANCE QUANTITATION

- 1. Extended Linear Dynamic Range**
- 2. Excellent Precision**
- 3. Excellent Accuracy**


# TSQ Quantum +ESI/SRM Quantitative Assay of Cabergoline

Used in the treatment of Parkinson's disease

Highly potent – typically present in extremely low levels in plasma





**Cabergoline**



**Pergolide  
(Internal Standard)**


# TSQ Quantum +ESI/SRM Quantitative Assay of Cabergoline in Bovine Plasma

5 $\mu$  Xterra MS18 150 x 2.1mm, ACN/NH<sub>4</sub>CH<sub>3</sub>COO (60/40, v/v), CE 19V, Q1 0.7 Da Q3 0.7 Da



# TSQ Quantum +ESI/SRM Quantitative Assay of Cabergoline in Bovine Plasma

150 x 2.1mm 5 $\mu$  Xterra MS18 LC Column, ACN/NH<sub>4</sub>CH<sub>3</sub>COO (60/40, v/v), SRM m/z 452 - 381  
Collision Energy 19V, Internal standard Pergolide, Q1 0.7 Da Q3 0.7 Da



**50 fg - 5 ng on column, 5 orders of Linear Dynamic Range - Excellent**

# TSQ Quantum

## +ESI/SRM Quantitative Assay of Cabergoline in Bovine Plasma


150 x 2.1mm 5 $\mu$  Xterra MS18 LC Column, ACN/NH<sub>4</sub>CH<sub>3</sub>COO (60/40, v/v), SRM m/z 452 - 381  
Collision Energy 19V, Internal standard Pergolide, Q1 0.7 Da Q3 0.7 Da

Nominal Amount (pg on column)	Mean Amount (pg on column)	Accuracy (% RE)	Precision (% CV)
0.050	0.057	13.6	9.3
0.500	0.508	1.6	8.4
5.000	4.709	-5.8	2.3
50.000	49.861	-0.3	1.0
500.000	462.354	-7.5	1.3
5000.000	5038.062	0.8	0.5

**Excellent Accuracy and Precision**


# TSQ Quantum, Low Femtogram +ESI/SRM at Unit Resolution

**Proprietary Pharmaceutical** in Mobile Phase, 3 $\mu$  Phenyl Hexyl 20 x 2mm,  
ACN/0.1%TFA/H<sub>2</sub>O/NH<sub>4</sub>CH<sub>3</sub>COO, +ESI, SRM m/z 425 - 218, Collision Energy 28V



# TSQ Quantum +ESI/SRM Quantitative Assay of Drug at Unit Resolution



Proprietary Drug in Mobile Phase, 3 $\mu$  Phenyl Hexyl 20 x 2mm,  
ACN/0.1%TFA/H<sub>2</sub>O/NH<sub>4</sub>CH<sub>3</sub>COO, SRM m/z 425 - 218, Collision Energy 28V



**1 fg - 1 ng on column (6 orders of linear dynamic range!)**


# High Performance Quantitation of a Drug in Bovine Plasma on the TSQ Quantum

Proprietary Drug, Unit Resolution, +ESI/SRM  $m/z$  644 - 114, Collision Energy 38V



# High Performance Quantitation of a Drug in Bovine Plasma on the TSQ Quantum

Proprietary Drug, Unit Resolution, +ESI/SRM m/z 644 - 114, Collision Energy 38V



# High Performance Quantitation of a Drug in Bovine Plasma on the TSQ Quantum

Proprietary Drug, Unit Resolution, +ESI/SRM m/z 644 - 114, Collision Energy 38V

## CALIBRATION STANDARDS - 8 replicates

Concentration	Area	ISTD Area	Area Ratio	Specified Amount (pg on column)	Calculated Amount (pg on column)	%Diff	%RSD
0.05pg/ $\mu$ L	14112	319466	0.0441	0.500	0.518	4%	3.3%
0.1pg/ $\mu$ L	25875	320259	0.0808	1.000	0.972	3%	3.3%
0.25pg/ $\mu$ L	61483	298452	0.206	2.50	2.52	-1%	2.3%
0.5pg/ $\mu$ L	128639	316553	0.406	5.00	5.00	0%	2.4%
1pg/ $\mu$ L	247948	313325	0.791	10.00	9.76	-2%	1.3%
2.5pg/ $\mu$ L	680210	335203	2.03	25.0	25.0	0%	1.3%
5pg/ $\mu$ L	1363255	332396	4.10	50.0	50.4	1%	1.1%
10pg/ $\mu$ L	2847221	347401	8.20	100.0	99.8	0%	1.2%

Excellent Precision and Accuracy

# High Performance Quantitation of a Drug in Bovine Plasma on the TSQ Quantum

Proprietary Drug, Unit Resolution, +ESI/SRM m/z 644 - 114, Collision Energy 38V


## Quality Control Samples - 8 replicates

Concentration	Area	ISTD Area	Area Ratio	Specified Amount (pg on column)	Calculated Amount (pg on column)	%Diff	%RSD
0.05pg/ $\mu$ L	13631	311163	0.04380	0.500	0.514	3%	3.4%
0.1pg/ $\mu$ L	26188	313428	0.08356	1.00	1.01	1%	2.8%
1pg/ $\mu$ L	251810	31593	0.79686	10.0	9.83	-2%	0.7%
10pg/ $\mu$ L	2649398	323837	8.18253	100	99.7	0%	2.2%

Excellent Accuracy and Precision Again

# High Performance Quantitation of a Drug in Bovine Plasma on the TSQ Quantum


Proprietary Drug, Unit Resolution, +ESI/SRM m/z 644 - 114, Collision Energy 38V



# TSQ Quantum

## Positive APCI/SRM of a Proprietary Drug in Bovine Plasma

$$Y = -0.00332178 + 0.0258183 * X \quad R^2 = 0.9993 \quad W: 1/X$$




**500fg – 5000pg on column ( $10^4$  dynamic range)**

# TSQ Quantum

## Negative ESI/SRM Quantitative Assay - $10^4$ Dynamic Range

Proprietary Drug, Bovine Plasma. Unit Resolution, SRM m/z 440 - 190, Collision Energy 32 V


$$Y = 308.156 + 3618.72 * X \quad R^2 = 0.9982 \quad W: 1/X$$



# TSQ Quantum

## Ruggedness Over 200 Injections in Plasma

**Peak area response of 200 injections of Paroxetine in protein precipitated plasma**



Thanks to Mark Churchill and Mark Harrison  
for data

**Thermo**Finnigan


# TSQ QUANTUM – HIGH RESOLUTION

**TSQ Quantum is the **Only**  
Triple Quadrupole Mass Spectrometer  
With High Resolution Capability**

**What is High Resolution?**


**How can High Resolution be Beneficial  
In Quantitation?**

# TSQ Quantum: Increasing Resolution For Reserpine (MW 608) in +ve ESI/MS




# TSQ Quantum: Increasing Resolution For Reserpine (MW 608) in +ve ESI/MS


**0.3 Da FWHM**



**0.2 Da FWHM**




**0.1 Da FWHM  
(High Resolution)**




**\*Signal Intensity at Unit Resolution = 2.1e6**

**Loss of intensity of ONLY x2-3 typically seen from Unit to High Resolution**


# Scanning Ion Stability for Quadrupoles



# TSQ Quantum - Instrument Tune Window Compound Optimization - Resolution Setting



# TSQ Quantum - Instrument Tune Window Compound Optimization - Resolution Setting



**High resolution achieved by a simple click of the mouse button**

# High Resolution on the TSQ Quantum

Quadrupole mass filtering fields are defined by:

1. electrode shape
2. field radius
3. frequency
4. voltage
5. length





The TSQ Quantum has true hyperbolic electrodes, large field radius, high frequency and voltage, and a long length (25 cm) and **hence has the ability to achieve high mass resolution**

# TSQ Quantum

## SRM Selectivity at High Resolution for Quantitation

Metabolite Quantitation, **Plasma extract**, MeOH/H<sub>2</sub>O, 500 ul/min,  
+ESI/SRM 372–328, Fast LC/MS





**Specificity of high resolution allows removal of two rogue peaks!!**

**Much less noise in high resolution**

# TSQ Quantum

## LC/ESI/SRM of Clenbuterol – Urine Matrix




Thanks to Mark Churchill and Mark Harrison

**Thermo**Finnigan

# TSQ Quantum - Unit vs. High Resolution +ESI/SRM Example of Great Sensitivity in HR

Pharmaceutical compound in Plasma, H<sub>2</sub>O/MeOH/ACN/0.1% formic,  
250 µl/min, 100 x 2 mm C8, +ESI, SRM m/z 311 - 259, Collision Energy 28V




**Low Resolution**

**10 femtograms  
on-column**

**Q1 0.7 Da FWHM  
Q3 0.7 Da FWHM**

**Peak Area 4144  
S/N 24**



**High Resolution**

**10 femtograms  
on-column**

**Q1 0.1 Da FWHM  
Q3 0.5 Da FWHM**


**Peak Area 1469  
S/N 24**

**Can use High Resolution Capabilities as low as 10 fg on-column!!**


# TSQ Quantum

## High Resolution Quantitative Assays

### +ESI/SRM Quantitation for a Metabolite of the Insecticide, Thiamethoxam, using High Resolution





**Thiamethoxam**  
MW 291



**Thiamethoxam Metabolite**  
MW 249

# +LC/ESI/SRM of Thiamethoxam Metabolite in Blueberry Extract TSQ Quantum, Unit versus High Resolution, 100 fg on column



C8 100 x 2.0 mm, ACN/0.1%FA:H<sub>2</sub>O/ 0.1%FA, SRM m/z 250-132, Collision Energy 14V



# TSQ Quantum: High Resolution +ESI/SRM Quantitation of Thiamethoxam Metabolite in Blueberry Extract

Blueberry Extract, C8 100 x 2.0 mm, Q1 0.2 Da Q3 0.5 Da, ACN/0.1%FA:H<sub>2</sub>O/ 0.1%FA, 100 fg – 1 ng on column - 4 orders of magnitude, SRM m/z 250-132, CE 14V

$$Y = 228.387 + 8365.93 * X \quad R^2 = 0.9994 \quad W: 1/X$$



**Lovely linearity**

# TSQ Quantum: High Resolution +ESI/SRM Quantitation of Thiamethoxam Metabolite in Blueberry Extract

**Matrix = Blueberry Extract, C8 100 x 2.0 mm, Q1 0.2 Da Q3 0.5 Da, ACN/0.1%FA:H<sub>2</sub>O/0.1%FA, SRM m/z 250-132, Collision Energy 14 V**

**Curve Index**   **Weighting Index**   **Origin Index**   **Equation**  
Linear                      1/X                      Ignore                       $Y = 228.387 + 8365.93 * X$    **R<sup>2</sup> = 0.9994**

<b>Amount</b>	<b>Area</b>	<b>Specified Amount (pg)</b>	<b>Calculated Amount (pg)</b>	<b>%RSD</b>
100fg on column	964	0.100	0.088	14.3%
1pg on column	9037	1.000	1.053	1.6%
10pg on column	85375	10.000	10.177	3.4%
100pg on column	882912	100.000	105.509	1.6%
1ng on column	8318237	1000.000	994.272	2.0%

# Conclusion

Because of the high resolution of the TSQ Quantum, we gain at least **an order in magnitude in quantitation limits**

Quantitation down to 100 fg of the metabolite on column is possible without interferences of the severe Blueberry matrix. This is **50 times lower** than the method proposed by the producer

The TSQ Quantum has over **4 orders of magnitude** a linear dynamic range in high resolution mode


# TSQ Quantum - Ruggedness

## Peak Area Stability at Unit (0.7 Da FWHM) and High Resolution (0.1 Da FWHM) with Time

### Clenbuterol Peak Area at 0.7 Da and 0.1 Da FWHM

### 250 injections over 45 hours

ThermoHypersil Keystone C18 100 x 2.1 mm  
Methanol / Ammonium acetate 5mM, 0.3 ml/min  
10  $\mu$ l injection of 10 pg/ $\mu$ l



Thanks to Mark Churchill and Mark Harrison  
for data

**Thermo**Finnigan

# TSQ Quantum – Qualitative Applications

## Tools For Structural Elucidation

Accurate Mass Measurement  
Clean High Quality MS/MS  
Different Scanning Functions


# TSQ Quantum – Qualitative Applications

## Tools For Structural Elucidation

**Accurate Mass Measurement**

# TSQ Quantum: Accurate Mass Determination of Parent Molecular Ion of Testosterone by LC/ESI/MS

## Testosterone, MW 288






PEG Internal Lock Masses, Mass Range scanned ~ 60 amu

# TSQ Quantum, Accurate Mass Determination of $MH^+$ , $m/z$ 289, of Testosterone by LC/ESI/MS (Q3 = 0.1 Da)


**Aver. Determined Mass = 289.2166, Proposed Exact Mass = 289.2168, Diff < 0.5 mmu**

**Internal PEG Lock Masses (\*)**



# TSQ Quantum: Accurate Mass Determination of Parent Molecular Ion of Erythromycin by LC/ESI/MS

## Erythromycin, MW 733





$\text{MH}^+ = \text{C}_{37}\text{H}_{68}\text{NO}_{13}$  Exact mass – 734.4688

PEG Internal Lock Masses, Mass Range scanned ~ 60 amu

# TSQ Quantum, Accurate Mass Determination of $MH^+$ , $m/z$ 734, of Erythromycin by LC/ESI/MS (Q3 = 0.1 Da)


## Internal PEG Lock Masses (\*)



**Aver. Determined Mass = 734.4696, Proposed Exact Mass = 734.4688, Diff < 1mmu**

# TSQ Quantum – Metabolite Identification through Accurate Mass Measurement

## 1-Hydroxymidazolam, Metabolite




Exact Mass  $\text{MH}^+ = 342.0809$


### Experiment

Accurate Mass Measurement of  $\text{MH}^+$  and Fragment Ion,  $[\text{MH}-18]^+$   
by LC/ESI/MS/MS using External PEG Lock Masses


# TSQ Quantum: Accurate Mass of $MH^+$ and $[MH-18]^+$ by LC/ESI/MS/MS using External PEG Lock Masses ( $m/z$ 305, 349) (6 LC inj, SIM, Q3 = 0.1 Da, Collision Gas = 1.2 mtorr, CE = 20 eV)



**First inj**  
**1**



**Last inj**  
**6**



**TSQ Quantum: Accurate Mass of MH<sup>+</sup> and [MH-18]<sup>+</sup> by LC/ESI/MS/MS  
using External PEG Lock Masses (m/z 305, 349)  
(6 LC inj, SIM, Q3 = 0.1 Da, Collision Gas = 1.2 mtorr, CE = 20 eV)**

	Time of Elution (min)	Accurate Mass MH <sup>+</sup> (u)	Standard Deviation (mmu)	Accurate Mass [MH-18] <sup>+</sup> (u)	Standard Deviation (mmu)
LC Inj 1	3.28	342.0827	0.3	324.0714	0.2
LC Inj 2	6.54	342.0832	0.3	324.0732	0.1
LC Inj 3	9.80	342.0830	0.1	324.0727	0.3
LC Inj 4	14.09	342.0828	0.3	324.0718	0.1
LC Inj 5	17.35	342.0820	0.4	324.0721	0.2
LC Inj 6	20.60	342.0817	0.1	324.0707	0.1
<b>Average</b>		<b>342.0826</b>		<b>324.0720</b>	
<b>Exact Mass</b>	<b>MH<sup>+</sup></b>	<b>342.0809</b>	<b>[MH-H<sub>2</sub>O]<sup>+</sup></b>	<b>324.0704</b>	


**Good agreement between exact masses of known MH<sup>+</sup>, proposed [MH-H<sub>2</sub>O]<sup>+</sup> and the average determined accurate masses (<2 mmu)**

# TSQ Quantum

## Metabolite Identification through Accurate Mass Measurement

**Typical situation - the metabolite (R) is unknown and metabolite identification is required through MS**

**However the active drug (M), in this case Midazolam, is known**



**Midazolam**  
**C<sub>18</sub>H<sub>13</sub>FN<sub>3</sub>Cl**  
**MW = 325**

**Since by ESI/MS it is known that the metabolite RH<sup>+</sup> is m/z 342 we know we have a [M+16] metabolite**

# TSQ Quantum

## Metabolite Identification through Accurate Mass Measurement

**Accurate Mass determined for the Metabolite,  $[M+16]H^+$   
on the TSQ Quantum = **342.0826****

Using an elemental composition calculator where the elemental lower limits are set by the composition of parent drug, Midazolam, (C = 18, H = 13, F = 1, N = 3, Cl = 1) the **ONLY** elemental compositions corresponding to the TSQ Quantum-determined accurate mass for  $[M+16]H^+$  (50 mmu window) are:

Exact mass of  $[M+O]H^+$  = 342.0809, **1.7 mmu different**

Exact  $[M+NH_2]H^+$  = 342.1048, **22.1 mmu different**

Exact  $[M+CH_4]H^+$  = 342.1173, **34.7 mmu different**

**Thus, Accurate Mass Measurement on the TSQ Quantum Provides  
Clear Identification of  $[M+16]H^+$  as the  $[M+O]H^+$  Metabolite**


# TSQ Quantum – Qualitative Applications

## Tools For Structural Elucidation

**Clean High Quality MS/MS**

# TSQ Quantum, Separation of Isobaric Compounds By High Resolution for Structural Elucidation

Mixture of Clenbuterol and 19-Norandrosterone, +ESI/MS, Q1 at Unit Resolution




Thanks to Mark Churchill and Mark Harrison  
for data

**Thermo**Finnigan

# TSQ Quantum, Separation of Isobaric Compounds By High Resolution for Structural Elucidation

Mixture of Clenbuterol and 19-Norandrosterone, +ESI/MS/MS, Q1 at Unit Resolution



Mixed product ion  
spectrum in unit resolution


Impossible for  
structural elucidation  
of unknowns

Thanks to Mark Churchill and Mark Harrison  
for data

**Thermo**Finnigan

# TSQ Quantum, Separation of Isobaric Compounds By High Resolution for Structural Elucidation

Mixture of Clenbuterol and 19-Norandrosterone +ESI/MS, Q1 at High Resolution




**Separation by High  
Resolution of Components  
only 0.13 Da Apart**

Thanks to Mark Churchill and Mark Harrison  
for data


**Thermo**Finnigan

# TSQ Quantum, +ESI/MS/MS, Q1 High Resolution Separation and MS/MS Characterization of Isobaric Species

Mixture of Clenbuterol and 19-Norandrosterone +ESI/MS/MS, Q1 at High Resolution



Clenbuterol  
MS/MS 277.1



19-norandrosterone  
MS/MS 277.2

Thanks to Mark Churchill and Mark Harrison  
for data


**Thermo**Finnigan

# TSQ Quantum – Qualitative Applications

## Tools For Structural Elucidation


### Different Scanning Functions

# Constant Neutral Loss Scanning



**Q3 is offset by the neutral loss under investigation**

**Only those compounds which give a fragment having that specific loss are detected**




# TSQ Quantum - Constant Neutral Loss Scanning

## Increases Specificity and Decreases Chemical Noise


Invitro Microsomal  
Incubation of Drug P  
(5  $\mu$ M)

Generate Metabolites with  
masses 16 and 34 da  
greater than parent drug



**Full-Scan MS**  
**50  $\rightarrow$  450 da**


No sign  
of metabolites




**Neutral Loss Scan**  
**50  $\rightarrow$  450 da**  
**NL = 116 da**

Clearly see  
metabolites

# Precursor Ion Scanning



# TSQ Quantum, Precursor Ion Analysis of Proprietary Drug and its Metabolites



# TSQ Quantum - Summary

# TSQ Quantum Project Team

## **Product Development**

Lance Arbini  
Paul Atherton  
Huy Bui  
Clay Campbell  
Pascual Cardenas  
Richard Currell  
Wayne Dewey  
Jean-Jacques Dunyach  
Phil Fong  
Jerry Glass  
Nigel Gore  
Trevor Hall  
Richard Hartford  
Rex Heller  
Steve Hurwitz  
Hugh Leonard  
Ray Lewis  
Jaime McGregor  
Teresa McLaughlin  
Iain Mylchreest

Terry Olney  
Carlos Rojas  
Alan Schoen  
Hans Schweingruber  
Bill Siebert  
Dennis Taylor  
Eugene Zhuk

## **Manufacturing**

Alvin Abdon  
Ken Bradbury  
Rey Brito  
Steve Brown  
Bob Cook  
Rudy Cruz  
Armando Delacruz  
Lenis Doan  
Gary Fontana  
Bret Johnson  
Joe Leon Guerrero  
Alex Manalo  
Frank Marcos

Misael Martinez  
Ed Matias  
Jose-Cruz Medina  
Gerry Rojas  
Phillip Tran

## **Marketing**

Gary Paul  
Mark Churchill  
Nelson Cooke  
Iain Green  
Mark Harrison  
Louis Maljers

## **Product Support**

Fred Ayres  
Ed Gonzales  
Bruce Hillestad  
**Technical Publications**  
Antony Harvey  
Rob Heather  
Toby Kraft  
Susan Lee  
Franklin Mason

## **Chemists**

Witold Winnik  
Mark Kagan  
Maurizio Splendore  
Liam Moran